

Evaluation of Small Mass Spectrometer Systems

C Richard Arkin¹, Timothy P. Griffin¹, Andrew K. Ottens²

Jorge A. Diaz³, Duke W. Follistein⁴

Fredrick W. Adams⁴, and William R. Helms⁴

¹ Dynacs Inc.

² University of Florida

³ Universidad de Costa Rica

⁴ NASA, Kennedy Space Center

Introduction

Understand Aspects of Designing a
Miniature Mass Spectrometer (MS) System

Systems Evaluated

Linear Quadrupole	Quadrupole Array
Quadrupole Ion Trap	Time of Flight
Sector	

Figures of Merit

Accuracy	Precision
Limits of Detection	Scan Rate
Response Time	Recovery Time
Volume	Weight

Set Up Scale to Rank Systems

Why Miniaturize?


- Improve Portability
Field Applications
Real-time Analysis
Decrease Contamination
- Size Reduction
Increase Sensor Density
Less Intrusive
- Cost Reduction
Less Weight
Less Power

Potential Applications

- Leak Detection
 - Hazardous Gases
 - Helium
 - Refrigerants
- Air Analysis
 - Worker Safety
 - Public Safety
 - Chemical Weapons
- Law Enforcement
 - Contraband Detection
- Field Medical Analysis
- Process Control
 - Semi-conductors
 - H₂ Economy


NASA Applications

- Next Generation Leak Detection
 - Leaks During Processing
 - Leaks During Cryogenic Fueling
 - Leaks During Launch Countdown
- Shuttle Engine Monitoring
- Air Analysis
 - International Space Station Air Lock
 - Shuttle Air Lock
- Process Control
 - Martian Fuel Generators



Current Problems


Response Time Scan Rate Sampling Density



Stanford Research Systems (SRS) XPR-2 Inficon XPR-2

RG-100

- | | |
|--------------------------------------|--------------------------------------|
| • Linear Quadrupole Analyzer | • Linear Quadrupole Analyzer |
| Cylindrical Rods: 6.35mm OD | Hyperbolic Rods |
| Rod Length: 11.4 cm | Rod Length: 12.7 mm |
| Inscribed Radius: xx | Inscribed Radius: 0.33 mm |
| Frequency: 2.76 MHz | Frequency: 13 MHz |
| • Pressure = 5×10^{-5} torr | • Pressure = 1×10^{-4} torr |





Ferran SRS & XPR Time-of-Flight (TOF) Ion Works

- Quadrupole Array System
16 Cylindrical Rods: 1 mm OD
Rod Length: 10 mm
 $\Omega = 16 \text{ MHz}$
- Reflectron TOF System
Orthogonal Acceleration
Flight Length: m


- Pressure: $\sim 5 \times 10^{-4}$ torr

- Pressure: $\sim 5 \times 10^{-6}$ torr




University of Florida UF-IT & P-Q Thermo Finnigan Ion Trap (UF-IT) Polaris-Q

- Quad. Ion Trap System
Ring Radius: 10 mm
Stretched Geometry
 $\Omega = 2.5$ MHz
No Buffer Gas
- Pressure: 4×10^{-6} torr




- Quad. Ion Trap System
Ring Radius: 7 mm
Stretched Geometry
 $\Omega = 1.03$ MHz
Helium Buffer Gas
- Pressure: $1 \times 10^{-3} (4 \times 10^{-6})$ torr



Monitor Group MG MG-2100

- Cycloidal Sector System
B: 0.5 Tesla
Pitch: 1 inch


- Pressure: $\sim 5 \times 10^{-6}$ torr



University of Minnesota CDFMS Compact Double Focus MS (CDFMS)

- Double Focus Sector System
B: 0.75 Tesla
Sector Radius: 20 mm


- Pressure: 10^{-5} torr




Experimental Parameters

Parameters	Definition	Shuttle Requirements
Accuracy	$\frac{[Test]_{meas} - [Test]_{true}}{[Test]_{true}} \times 100\%$	< 10% or 5ppm, whichever is greater
Precision	$\frac{\sigma_{Test, meas}}{[Test]_{true}} \times 100\%$	< 5% or 3ppm, whichever is greater
Limit of Detection	$3 \sigma_{zero}^a$	H ₂ , O ₂ : 25 ppm; He: 100 ppm; Ar: 10 ppm ^b
Response Time	Time required for response from valve change to 95% new reading	10 s
Recovery Time	Time after valve change to reach 5% of previous sample reading (new sample is zero)	30 s
Scan Time	Experiment Time / Scans	1 s
Maximum Size	Sum of Individual Components	$3.5 \times 10^4 \text{ cm}^3$
Maximum Weight	Sum of Individual Components	10 kg
^a Theoretical limit of detection		
^b Measured limit of detection		


Experimental Scan




Accuracy




Precision




Limits of Detection (LOD)




Response Time




Recovery Time




Scan Time



System Volume



System Weight



Evaluation Chart

	Accuracy	Precision	LOD	Response	Recovery	Scan Rate	System Volume	System Weight	Score
SRS	2	2	1	2	6	8	6	7	4.3
XPR-2	3	8	2	7	5	5	2	6	4.8
Ferran	10	10	7	10	10	9	3	6	8.1
Polaris-Q	7	3	2	8	8	2	10	10	6.3
UF-IT	3	3	6	3	3	2	7	8	4.4
TOF	9	6	7	9	7	2	9	8	7.1
MG-2100	9	7	10	10	5	10	6	7	8.0
CDFMS	4	5	9	8	6	2	2	4	5.0
Average ¹	4.7	4.5	4.5	6.2	5.8	3.5	6.0	5.5	

¹Excluding Ferran & MG-2100

Conclusions

- Various Mass Analyzer Systems Evaluated
- Several Systems Show Promise
 - Stanford Research Systems RGA-100
 - Inficon XPR-2
 - University of Florida – Ion Trap
 - Compact Double Focus Mass Spectrometer
- Areas that Need Improvement
 - Response Time Recovery Time
 - System Volume System Weight
- Future Work
 - Investigate Techniques to Improve Systems
 - Evaluate Engineering Challenges